

## Product Specifications

<b>Product number:</b>	<i>BSE-1000</i>
<b>Description:</b>	<i>Bovine spleen extract, activities include Ro (SS-A), La (SS-B), SCL-70, Sm, and RNP.</i>
<b>Amount:</b>	<i>1 ml</i>
<b>Physical state:</b>	<i>Lyophilized powder, buffered with 0.01 M phosphate, pH 7.2, containing 0.02% NaN<sub>3</sub>.</i>
<b>Manufacturing method:</b>	<i>Saline soluble material is extracted from bovine spleens with physiological buffer. Insoluble material is removed by centrifugation and filtration. The extract is optimized for autoantigen activity in immunodiffusion assays.</i>
<b>Storage conditions:</b>	<i>&lt;-40°C</i>
<b>Instructions for use:</b>	<i>Reconstitute using 1 ml of deionized water. Allow several minutes for dissolution. Vortexing will speed dissolution.</i>
<b>Warnings:</b>	<i>For research or further manufacturing use only.</i>