

## **Product Specifications**

<b>Product number:</b>	<i>ENA-1001</i>
<b>Description:</b>	<i>ENA Extract. Precipitates antisera to SSA (Ro), SSB (La), Scl-70, Sm, RNP, and Jo-1.</i>
<b>Amount:</b>	<i>0.25 ml</i>
<b>Physical state:</b>	<i>Lyophilized. Buffered with 0.01M phosphate in saline.</i>
<b>Manufacturing method:</b>	<i>Bovine and rabbit tissues are extracted in phosphate buffer. Insoluble material is removed. Extractable nuclear antigens are lyophilized.</i>
<b>Storage conditions:</b>	<i>Store <math>\leq 0^{\circ}</math> C. Aliquot and store reconstituted extract, do not subject to repeated freeze thaw.</i>
<b>Instructions for use:</b>	<i>Reconstitute with 0.25 ml cold distilled water. Allow 1 hr for complete dissolution. Vortexing may be required. Some assay systems may require more concentrated or more dilute extract. Some precipitate may occur, but it will not affect the product.</i>
<b>Warnings:</b>	<i>For research or further manufacturing use only.</i>